

PRODUCT INFORMATION

Item Number	PIP25
Volume	350ml
Dimensions	75mm (H) x 157mm / 104mm (W)

FEATURES

Maceratable: Tested in accordance with PAS29:1999

Made from sustainable biodegradable cellulose materials
NHSSC NPC: FBF85178

Reduced cradle to grave carbon footprint compared to collection with single use plastic items

Reduced use of sterile single use plastic items

The high temperature manufacturing process demonstrates equivalency in urine MSU culture results compared to sterile single use plastic items.

Product stackable and easy to separate, reducing storage space.

Single use thermofibre urine collection device, for transferring to the correct primary tube for laboratory MSU analysis.

SECTORS

Medical • Pharmaceutical and Laboratory

STORAGE INSTRUCTIONS

Store in dry conditions (not a sluice room), away from direct sunlight and extreme heat. For ward use. You are advised to retain the packaging for reference. Always inspect product prior to use and do not use if damaged.

PACKING

25 pots per box
10 boxes per case



INNER PACKAGING BARCODE

5024951030928

OUTER PACKAGING BARCODE

5024951031147

For further technical information, please contact Customer Services on 0333 320 8550 or the Technical team on 0333 320 8555